

TABLE 1: AIRWAY DEVICES

Age (years)	Weight (kg)	Guedel		LMA	iGel		Laryngoscope Blade/CMAC	Air-Q (max size uncuffed ETT)
Neonate	1	3.5	5	1	1		Mac 0 / Miller 0	0 (3) 0.5 (4)
	2							
	3							
	4							
Infant	5	5.5		1.5	1.5		Mac 1 / Miller 1	1.0 (4.5) 1.5 (5)
	6							
	7							
	8							
1 to 6	9	6.5 (for 1-6yr old)		2	2		Mac 2	2 (5.5)
	10							
	11							
	12							
	13							
	14							
	15							
	17							
6 to 10	20	7		2.5	2.5		Mac 2 or 3	2 (5.5)
	25							
	30							
10 to 16	35	8 (12yrs to adult)		3	3		Mac 3	3 (7)
	40							
	45							
	50							
	55							
	60							

TABLE 2: ET TUBES: <1y SIZE ON WEIGHT, >1Y SIZE ON AGE

Age (years)	Weight (kg)	ETT Size (Microcuff)	ETT Size (Uncuffed)	Oral ETT Length (cm)	
Neonate	1		2 - 2.5	7	
	2		2.5 - 3	8	
	3		3 - 3.5	9	
	4		3.5	10	
Infant	5	3.5		11	
	6		3.5-4	11.5	
	7	12			
	8	12.5			
	9	13			
1		3.5	4	13.5	
2			4	4.5	14
3				4.5-5	14.5
4		4.5	5	15	
5			5-5.5	15.5	
6		5	5.5	16	
7			5.5-6	16.5	
8		5.5	6	17	
9				17.5	
10		6		18	
11				18.5	
12		6.5		19	
13				19.5	
14		7		20	
15				20.5	
ALWAYS HAVE ONE LARGER AND ONE SMALLER AVAILABLE					
* Microcuff ETT should only be used in term neonates >3kg					

ALWAYS HAVE ONE LARGER AND ONE SMALLER AVAILABLE
*** Microcuff ETT should only be used in term neonates >3kg**
Please note: No bars on the LMA to allow passage of a tube through it. Always check that a tube can be passed first.

TABLE 3: DIFFICULT AIRWAY/ACCESSORY DEVICES: SIZE BASED ON ETT INTERNAL DIAMETER

ETT Size ID	Stylet	Bougie	Suction Catheter	NGT	Fibreoptic Scope		Video Bronchoscope	Cook Airway Exchange Catheter											
3	Small	5 Fr	6 Fr	6 Fr		Olympus 2.2mm	Size 3 (No 162) *SNUG, Remove 15mm Connector*	8 Fr											
3.5		10 Fr			8 Fr		8 Fr		Storz 2.8mm	Olympus 2.2mm	Size 3 (No. 162)								
4			8 Fr	8 Fr		Storz 2.8mm		Olympus 2.2mm				Size 3 (No. 162)							
4.5					Medium		15 Fr		10 Fr	10 Fr	Storz 4.0mm		Olympus 4.1mm	Size 4 (No. 164)					
5	10 Fr	10 Fr	Storz 4.0mm	Olympus 4.1mm		Size 4 (No. 164)													
5.5								12 Fr	12 Fr	Storz 4.0mm		Olympus 4.1mm			Size 4 (No. 164)				
6	12 Fr	12 Fr														Storz 4.0mm	Olympus 4.1mm	Size 4 (No. 164)	
6.5								14 Fr	14 Fr										Storz 4.0mm
7	14 Fr	14 Fr																	
7.5					14 Fr		14 Fr	Storz 4.0mm	Olympus 4.1mm		Size 4 (No. 164)								

TABLE 4: INTUBATING THROUGH 1st GENERATION DISPOSABLE LMA

LMA	Uncuffed ETT	Cuffed ETT
1	3	LMA under size 3 will not allow passage of pilot balloon
1.5	3.5	
2	4	
2.5	5	
3	6	
		5.5